EN	ORTHOFIX® GALAXY FIXATION SYSTEM 2	JA	オーソフィックス® GALAXY 固定システム 111
IT	ORTHOFIX® GALAXY FIXATION SYSTEM 11	ZH	ORTHOFIX® GALAXY 固定支架系统 122
FR	SYSTÈME DE FIXATION GALAXY D'ORTHOFIX® 20	CS	FIXAČNÍ SYSTÉM ORTHOFIX® GALAXY 131
DE	ORTHOFIX® GALAXY-FIXATIONSSYSTEM 29	PL	ORTHOFIX® SYSTEM STABILIZACJI GALAXY 140
ES	ORTHOFIX® SISTEMA DE FIJACIÓN GALAXY 38	SL	ORTHOFIX® PRITRDILNI SISTEM GALAXY 149
DA	ORTHOFIX® GALAXY-FIXATIONSSYSTEM 47	KO	ORTHOFIX® GALAXY 고정 시스템 158
FI	ORTHOFIX® GALAXY -KIINNITYSJÄRJESTELMÄ 56	AR	ORTHOFIX® GALAXY FIXATION SYSTEM 168
NO	ORTHOFIX® GALAXY-FIKSERINGSSYSTEM 65	BG	СИСТЕМА ЗА ФИКСАЦИЯ GALAXY HA ORTHOFIX® 177
NL	ORTHOFIX® GALAXY FIXATIESYSTEEM 74	LT	„ORTHOFIX®“ FIKSAVIMO SISTEMA „GALAXY“ 187
PT	ORTHOFIX® SISTEMA DE FIXAÇÃO GALAXY 83	MS	ORTHOFIX® SISTEM FIKSASI GALAKSI 196
SV	ORTHOFIX® GALAXY FIXERINGSSYSTEM 92	RO	SISTEMUL DE FIXARE ORTHOFIX® GALAXY 205
EL	ORTHOFIX® ΣΥΣΤΗΜΑ ΣΤΑΘΕΡΟΠΟΙΗΣΗΣ GALAXY 101	TR	ORTHOFIX® GALAXY FİKSASYON SİSTEMİ 214

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Telephone 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

PQ GAL L 06/19

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online.

EN

Important information - please read prior to use

ORTHOFIX® GALAXY FIXATION SYSTEM

Orthofix Srl
Via delle Nazioni 9 - 37012 Bussolengo (VR) Italy
 Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

PRODUCT INFORMATION

DESCRIPTION

The Orthofix Galaxy Fixation System consists of a series of components used in conjunction with Orthofix bone gripping elements (bone screws, wires...). External fixation systems are modular, therefore different frame configurations are possible.

The Orthofix Galaxy Fixation System components are not intended to replace normal healthy bone or to withstand the stresses of full weightbearing, particularly in unstable fractures or in the presence of non union, delayed union or incomplete healing. The use of external supports (e.g. walking aids) is recommended as a part of the treatment. The system consists of various modules to be applied in different anatomical sites of the upper and lower limb as well as the pelvis. When used correctly, the Orthofix Galaxy Fixation System maintains limb function, minimizes surgical trauma to anatomical structures, preserves the blood supply and osteogenic potential of the tissues. All Orthofix devices are intended for professional use only. Surgeons who supervise the use of Orthofix devices must have full awareness of orthopaedic fixation procedures as well as adequate understanding of the philosophy of the Orthofix modular system. To promote the proper use of its fixation system, and establish an effective promotional and training tool, Orthofix has developed several manuals or CD-ROM's containing the relevant information (i. e. general philosophy, surgical application, etc.) called "Operative Techniques". These are available in several languages as a complimentary service for surgeons who have adopted the Orthofix system. If you wish to receive a personal copy, please contact Orthofix or its local authorized representative, with a description of the medical device to be used.

MATERIALS

The Orthofix Galaxy Fixation System is comprised of stainless steel, aluminium alloy, titanium alloy and plastic components. Those components which contact the patient are the percutaneous pins (bone screws), threaded wires, drill bits, guides used during screw insertion. These are manufactured from surgical grade stainless steel. Some of the Orthofix external fixation bone screws (pins) are supplied with a thin, plasma sprayed coating of hydroxyapatite (HA) on the threaded portion of the shaft.

INDICATIONS FOR USE

The Galaxy Fixation System is intended to be used for bone stabilization in trauma and orthopedic procedures, both on adults and all pediatric subgroups excepts newborns as required.

The indications for use include:

- open or closed fractures of the long bones.
- vertically stable pelvic fractures or as a treatment adjunct for vertically unstable pelvic fractures.
- infected and aseptic non-unions.
- joint pathologies/injuries of upper and lower limb, such as:

- proximal humeral fractures
- intra-articular knee, ankle and wrist fractures
- delayed treatment of dislocated and stiff elbows
- chronic, persistent elbow joint instability
- acute elbow joint instability after complex ligament injuries
- unstable elbow fractures
- additional elbow stabilization of post-operative unstable internal fixation

The Orthofix Galaxy Wrist external fixator is intended for the following indications:

- intra-articular or extra-articular fractures and dislocations of the wrist with or without soft tissue damage
- polytrauma
- carpal dislocations
- unreduced fractures following conservative treatment
- bone-loss or other reconstructive procedures
- infection

NOTE: The Shoulder Fixation System is intended to be used for proximal humeral fractures where two thirds of the metaphysis is intact.

CONTRAINdications

The Orthofix Galaxy Fixation System is not designed or sold for any use except as indicated.

Use of the system is contraindicated in the following situations:

- Patients with mental or physiological conditions who are unwilling or incapable of following postoperative care instructions.
- Patients with severe osteoporosis*.
- Patients with severe, poorly controlled diabetes mellitus.
- Patients with compromised vascularity.
- Patients with previous infections.
- Patients with malignancy in the fracture area.
- Patients with neuromuscular deficit or any other conditions that could influence the healing process.
- Patients who are HIV positive.
- Patients with foreign body sensitivity. Where material sensitivity is suspected, tests should be made prior to implant insertion.

* As defined by the World Health Organization: "Bone mineral density of 2.5 standard deviations or more below the mean peak bone mass (average of young, healthy adults) in presence of one or more fragility fractures".

WARNINGS & PRECAUTIONS

1. Fracture stabilization has to be done following correct fracture reduction.
2. The clamp must be closed first manually by turning the metal ring clockwise before locking it firmly by tightening the cam with the Universal T-Wrench or the 5mm Allen Wrench.
3. Large Clamp (93010), Medium Clamp (93110), Small Clamp (93310), Wrist Module (93350), Wire Locking Clamp (93620), Elbow Hinge (93410), Large-Medium Transition Clamp (93030), Small Multiscrew Clamp-Long (93320), Small Multiscrew Clamp-Short (93330), cannot be disassembled.
4. Before fixator application, make sure that the clamps are fully loosened.
5. Frame stability must be checked intra-operatively before the patient leaves the operating theatre.
6. Compression is never recommended in a fresh fracture.
7. Particular care should be taken to avoid that screws and wires enter the joints or damage the growth plates in children.
8. Depending on the clinical and radiological findings, the surgeon will decide on the number of rods and bone screws needed to achieve the appropriate frame stability.
9. Any device implanted into the patient, such as bone screws and threaded wires, and in general any device which is labelled "single use only": MUST NOT BE RE-US ED
10. Screw length and thread length should be selected in accordance with bone and soft tissue dimensions. Excessive penetration of the second cortex by any type of screw should be avoided, because of the risk of soft tissue damage.
11. The maximum diameter of the screw thread should not be greater than one third of the bone diameter (e.g. 6-5 or 6-5.6mm bone screws for bone diameter greater than 20mm).
12. For pre-drilled bone screws, pre-drilling with appropriate drill bits and drill guides prior to screw insertion is imperative. Matching grooves on screws and drill bits help the surgeon to use the correct drill bit. Blunt drill bits can cause thermal damage to the bone and should always be discarded.
13. Self-drilling screws with a thread diameter of 5.00mm or above should never be inserted with a power tool, but always by hand or with a hand drill. Self-drilling screws with smaller thread diameters may be inserted with a power drill at low speed.
14. Transfix pins of 6mm in diameter are self-drilling and may be inserted with a power drill. These pins are used in association with the Fixator for temporary ligamentotaxis of the ankle and knee. Orthofix Transfix pins are single use devices and should never be re-used. They are connected to the Galaxy Rods with two Large Clamps.
15. For more stable fixation of a fracture with a fixator, we recommend that the nearest bone screw is applied fairly close to the fracture margin (a minimum of 2cm is recommended) and that these distances are equal on both sides of the fracture.
16. Do not introduce non "MR"-marked devices in the MR environment.
17. Appropriate Orthofix instrumentation should be used to insert bone screws correctly.
18. All equipment should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
19. The fixator should be applied at a sufficient distance from the skin to allow for post-operative swelling and for cleaning, remembering that the stability of the system depends upon the bone-fixator distance. If the fixator is sited at a distance of more than 4cm from the bone the surgeon will decide on the number of rods and bone screws needed to achieve the appropriate frame stability.
20. Components may not be interchangeable between all Orthofix Fixation Systems. Consult individual operative technique guides for interchangeable components.
21. Additional equipment may be required for fixation application and removal such as screw cutters and power drill.
22. Screw and frame integrity should be monitored at regular intervals. To avoid causing injury, the end of threaded wires, bone screws and transfix pins that have been cut, should be protected with special covers.
23. Meticulous screw site hygiene is required.
24. All patients must receive instruction on the use and maintenance of their external fixation frame, and on pin site care.
25. Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
26. The fracture site gap should be reassessed periodically during healing, and adjustments to the frame made as necessary. Persistent separation of the fracture ends may lead to delay in bone union.
27. Removal of the device: the surgeon should make the final decision whether a fixation device can be removed.
28. Do not use components of the Orthofix Galaxy Fixation System in conjunction with products of other manufacturers, unless otherwise specified, as the combination is not covered by the necessary validation.

29. Elbow Distractor

The Elbow Distractor: is intended to be used to distract the joint intra-operatively in case of elbow stiffness.

- The entity of distraction of the Elbow must be verified under image intensification.
- It is mandatory to expose the Ulnar nerve prior to distraction.

30. Shoulder Fixation System

The tip of the threaded wire should be in the subchondral area of the humeral head.

- During wire insertion, use the wire guide to avoid damage to the wire and soft tissues and/or joint impingement. After wire insertion, check joint function.
- Insert Wires in the safe corridors to avoid damage to anatomical structures.
- Do not drill Wires into the soft tissues but push them through the skin. Use the drill at low speed when inserting the wires into the bone.
- 2.5mm threaded wires are used with the blocking wire clamp.
- The first threaded wire should be inserted in the centre of the humeral head in order to target its apex.
- Wires are cylindrical and can be backed out if necessary.

- Use the dedicated Orthofix instrumentation to insert threaded wires.
- The skin around the threaded wires (wires sites) must be meticulously cleaned.

31. Galaxy Wrist

- To ensure the correct locking of the multi-screw clamp use always 2 screws and ensure that they are of the same diameter.
- Depending on the fracture patterns additional fixation techniques may be necessary.

POSSIBLE ADVERSE EFFECTS

1. Nerve or vessel damage resulting from insertion of wires and screws.
2. Superficial or deep bone screw tract infection, osteomyelitis, or septic arthritis, including chronic drainage of bone screw sites after device removal.
3. Oedema or swelling; possible compartment syndrome.
4. Joint contracture, subluxation, dislocation or loss of range of motion.
5. Premature bone consolidation during distraction osteogenesis.
6. Failure of bone to regenerate satisfactorily, development of non union or pseudarthrosis.
7. Fracture of regenerate bone or through bone screw holes after device removal.
8. Loosening or breakage of implants.
9. Bony damage due to inappropriate implants selection.
10. Bone deformity or equinus of the foot.
11. Persistence or recurrence of the initial condition requiring treatment.
12. Reoperation to replace a component or entire frame configuration.
13. Abnormal growth plate development in patients who are skeletally immature.
14. Foreign body reaction to implants or frame components.
15. Tissue necrosis secondary to implants insertion.
16. Pressure on the skin caused by external components when clearance is inadequate.
17. Limb length discrepancy.
18. Excessive operative bleeding.
19. Intrinsic risks associated with anesthesia.
20. Intractable pain.
21. Bone sequestration secondary to rapid drilling of bony cortex with heat build-up and bone necrosis.
22. Vascular disorders including thrombophlebitis, pulmonary embolus, wound hematomas, avascular necrosis.

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure which require further surgical intervention to remove or replace the external fixation device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the external fixation devices are important considerations in the successful utilization of Orthofix external fixation devices by the surgeon. Proper patient selection and the patient's ability to comply with physician instructions and follow prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindications or is predisposed to any contraindications, DO NOT USE Orthofix Galaxy Fixation devices.

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

Implantable Device*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol "⊗" reported on the product label. After the removal from the patient, the implantable device* has to be discarded.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

Non Implantable Device

The "SINGLE USE" non implantable device of Orthofix is identified through symbol "⊗" reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions were written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the health care facility to ensure that the reprocessing is performed in accordance with instructions and to perform validation and routine monitoring of the process. Any deviation from these instructions is responsibility of the health care facility in charge of reprocessing.

STERILE & NON-STERILE Product

Orthofix provides certain devices STERILE while others are provided NON-STERILE.

Please review the product label to determine the sterility of each device.

Sterile

Devices provided STERILE are labeled as such.

Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.

Non-Sterile

Unless otherwise noted, Orthofix devices are provided NON-STERILE. Orthofix recommends that all NON-STERILE devices be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Product integrity and performance are assured only if packaging is undamaged.

Cleaning, Sterilization and Maintenance

Cleaning is an essential pre-requisite to ensure effective disinfection or sterilization. All instruments must be cleaned prior to use and after each re-use of the instruments. Cleaning may be done following the validated cleaning process described below and using washing instruments and machines validated by the hospital.

WARNINGS

- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning solutions with a pH 7-9.5 are recommended. Cleaning solutions with an higher pH value must be avoided.
- Aluminium based instruments are damaged by alkaline (pH>7) detergents and solutions.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used. The contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing in order to remove soiling that accumulates in recesses.

Limitations and restrictions on reprocessing

- Repeated reprocessing has minimal effect on reusable fixators and instruments.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only must not be reused regardless any reprocessing.

Point of use

- It is recommended that instruments are reprocessed as soon as is reasonably practicable following use, putting them inside a sink filled with cold water (< 40°C) or a neutral pH solution for at least 10 minutes and removing gross soiling with a soft cloth or a soft brush.
- Don't use a fixating detergent or hot water as this can cause the fixation of residue which may influence the result of the reprocessing process.

Containment and transportation

- Follow hospital protocols when handling contaminated and bio-hazardous materials. Used instruments should be covered to minimize the risk of cross contamination.
- All used surgical instruments should be regarded as contaminated. Their handling, collection and transportation must be strictly controlled to minimize any possible risks to patients, personnel and any area of the healthcare facility.

Preparation for decontamination

- Disassemble the devices where applicable. See Orthofix operative techniques for further details.

Cleaning: manual

- 1) Soak the single components in the cleaning solution. Orthofix recommends the use of a pH-neutral enzymatic cleaning agent. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- 2) Thoroughly scrub the single components in the cleaning solution with a soft brush until all visible soiling is removed. Metal brushes should be avoided. Use a soft brush to remove residuals from lumens using a twisting motion and use a syringe filled with detergent if necessary. When cleaning hinges make sure that all areas are reached.
- 3) Rinse the single components in running tap water.
- 4) Ultrasonic clean the single components in a degassed cleaning solution. Orthofix recommends the use of a pH-neutral enzymatic or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- 5) Rinse the components in sterile or freshly prepared purified water.
- 6) Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer for at least five minutes.

Cleaning: Automated

When the devices to be cleaned have lumens or present complexity, a preliminary manual cleaning may be required.

Use a validated, properly maintained and calibrated washer disinfector.

- 1) Place all the instruments into washer baskets:
 1. Place heavier devices in the bottom of the baskets
 2. Connect cannulations to the proper injector jets
 3. Wherever possible, all parts of disassembled devices should be kept together in one container
- 2) Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
- 3) Orthofix recommends the use of a pH-neutral enzymatic or slightly alkaline solution. When using alkaline solutions, a neutralizer must be added. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- Orthofix recommends that cycle steps are at least as follows:
 1. Prewash at low temperatures
 2. Main wash at 40-60°C for at least 5 minutes
 3. Rinse with demineralized water
 4. Thermal disinfection at 90-95°C for at least 5 minutes
- 4) Rinse the components in sterile or freshly prepared purified water.
- 5) Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Maintenance, inspection and testing

All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed.

All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.

Particular attention should be given to:

- Cannulated devices (NB: cannulated drill-bits are single-patient use only).
 - Cutting edges: Discard blunt or damaged instruments.
 - Hinged instruments: check for smooth and regular movement of hinges
 - Locking mechanisms should be checked for action.
 - Cutting instruments must be checked for sharpness.
 - When instruments form part of an assembly, check assembly with matching components.
 - Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil.
- The cams and bushes in the ball joints of the Elbow Distractor must be replaced after every use.

Packaging:

- Wrap the tray before sterilization with an approved sterilization wrap or insert it into a rigid sterilization container in order to prevent contamination after sterilization.
- Do not include additional systems or instruments in the sterilization tray. Sterility cannot be guaranteed if the sterilization tray is overloaded.
- The total weight of a wrapped instrument tray should not exceed 10kg.

Sterilization:

- Steam sterilization is recommended. Gas plasma, dry heat and EtO sterilization must be avoided as they were not validated for Orthofix products.
- Use a validated, properly maintained and calibrated steam sterilizer.
- The steam quality must be appropriate for the process to be effective.
- Do not exceed 140°C (284°F).
- Do not stack trays during sterilization.
- Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

STEAM STERILIZER TYPE	GRAVITY	PRE-VACUUM	PRE-VACUUM (Not recommended for use in US)
Minimum exposure Temperature	132°C (270°F)	132°C (270°F)	134°C (273°F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes
Drying Time	30 minutes	30 minutes	30 minutes

Storage

Store the sterilized instrument in the sterilization packaging in a dry and clean-environment at room temperature.

Disclaimer: "The instructions provided above have been validated by Orthofix as being a true description of the preparation of a device for first clinical use or for re-use of multiple use devices. It remains the responsibility of the reprocessor to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes should be adequately recorded. Likewise any deviation by the reprocessor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences and should also be appropriately recorded".

Note: ANY DEVICE WHICH IS LABELLED "SINGLE USE ONLY" MUST NEVER BE REUSED. ORTHOFIX IS ONLY RESPONSIBLE FOR SAFETY AND EFFECTIVENESS FOR THE FIRST PATIENT USE OF SINGLE USE DEVICES. The institution or practitioner bears full responsibility for any subsequent use of these devices.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Manufacturer contact:

Please contact your local Orthofix sales representative for further details.

MRI SAFETY INFORMATION

GALAXY WRIST

Non-clinical testing has demonstrated that the Galaxy Wrist Components are MR Conditional. A patient with the Galaxy Wrist components can be safely scanned under the following conditions:

- Static magnetic field of 1.5Tesla and 3.0Tesla.
- Maximum spatial magnetic field gradient of 900-Gauss/cm (90mT/cm).
- Maximum MR System reported-whole-body-averaged specific absorption rate (SAR) of <4.0W/kg (First Level Controlled Operating Mode).
- No local transmit/receive coils must be used on the device.
- The Galaxy Fixation System must be entirely outside the MR scanner bore. No part of the Galaxy Fixation System must extend into the MR bore. Therefore MR scanning of body parts where the Galaxy Fixation System is located is Contraindicated.

Heating Information

Under the scan conditions defined above, the Galaxy Fixation System frames are expected to produce a maximum temperature rise of 1°C after 15 minutes of continuous scanning.

Displacement Information

The system will not present an additional risk or hazard to a patient in the 1.5Tesla and 3Tesla MR environment with regard to translational attraction or migration and torque.

GALAXY FIXATION SYSTEM

Galaxy System Fixator Components are labeled MR CONDITIONAL  according to the terminology specified in ASTM F2503 Standard Practice for Marking Medical Devices and Other Items in the Magnetic Resonance Environment.

Non-clinical testing has demonstrated that the Galaxy System Fixator Components is MR Conditional according to the terminology specified in ASTM F2503 Standard Practice for Marking Medical Devices and Other Items in the Magnetic Resonance Environment. Non-clinical testing, done according to ASTM F2052-06, F2213-06, F2182-11, F2119-07, demonstrated that a patient with the Galaxy Fixation System can be safely scanned under the following conditions:

- Static magnetic field of 1.5Tesla and 3.0Tesla.
- Maximum spatial magnetic field gradient of 900-Gauss/cm (90mT/cm).
- Maximum whole-body-averaged specific absorption rate (SAR) of 4.0W/kg in the First Level Controlled Mode for 15 minutes of scanning.
- No local transmit/receive coils must be used on the device.
- The Galaxy Fixation System must be entirely outside the MR scanner bore.

No part of the Galaxy Fixation System must extend into the MR bore.

Therefore MR scanning of body parts where the Galaxy Fixation System is located is Contraindicated.

Displacement Information

The system will not present an additional risk or hazard to a patient in the 1.5Tesla and 3Tesla MR environment with regard to translational attraction or migration and torque.

Heating Information

Comprehensive electromagnetic computer modeling and experimental testing was performed on the following systems:

- 1.5-Tesla/64-MHz: Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner.
- 3-Tesla/128-MHz: Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI, Active-shielded, horizontal field scanner.

to determine the worst heating in seven configurations of Orthofix Galaxy Fixation System. From these studies, it is concluded that once the entire external fixation frame is visible outside the MRI bore, the maximum heating is less than 2°C. In non-clinical testing the worst scenarios produced the following temperature rises during MRI under the conditions reported above:

	1.5Tesla System	3.0Tesla System
Galaxy Fixation System		
Minutes of scanning	15	15
Calorimetry measured values, whole body averaged SAR (W/kg)	2.2W/Kg	2.5W/Kg
Highest temperature Rise less than (°C)	2°C	2°C

Please note that temperature changes reported apply to the designed MR systems and characteristics used. If a different MR system is used, temperature changes may vary but are expected to be low enough for safe scanning as long as all Galaxy System Fixator Components are placed **outside** the MR bore.

MR PATIENT SAFETY

MRI in patients with Galaxy Fixation System can only be performed under these parameters. It is not allowed to scan the Galaxy Fixation System directly. Using other parameters, MRI could result in serious injury to the patient. When the Galaxy Fixation System is used in conjunction with other External Fixation Systems please be advised that this combination has not been tested in the MR environment and therefore higher heating and serious injury to the patient may occur. Because higher *in vivo* heating cannot be excluded, close patient monitoring and communication with the patient during the scan is required. Immediately abort the scan if the patient reports burning sensation or pain.

Galaxy Fixation System can only be guaranteed for MRI when using the following components to build a frame:

(*the following components are listed in non-sterile configuration. Please consider that the same MRI information and performance are applicable to the same components in gamma-sterile configuration if available (code number preceded by 99-, e.g. 99-93030))

RODS*

Code	Description
932100	Rod 100mm long, 12mm diameter
932150	Rod 150mm long, 12mm diameter
932200	Rod 200mm long, 12mm diameter
932250	Rod 250mm long, 12mm diameter
932300	Rod 300mm long, 12mm diameter
932350	Rod 350mm long, 12mm diameter
932400	Rod 400mm long, 12mm diameter
99-932450	Rod 450mm long, 12mm diameter, sterile
99-932500	Rod 500mm long, 12mm diameter, sterile
99-932550	Rod 550mm long, 12mm diameter, sterile
99-932600	Rod 600mm long, 12mm diameter, sterile
99-932650	Rod 650mm long, 12mm diameter, sterile
939100	Rod 100mm long, 9mm diameter
939150	Rod 150mm long, 9mm diameter
939200	Rod 200mm long, 9mm diameter
939250	Rod 250mm long, 9mm diameter
939300	Rod 300mm long, 9mm diameter
936060	Rod 60mm long, 6mm diameter
936080	Rod 80mm long, 6mm diameter
936100	Rod 100mm long, 6mm diameter
936120	Rod 120mm long, 6mm diameter
936140	Rod 140mm long, 6mm diameter
936160	Rod 160mm long, 6mm diameter
936180	Rod 180mm long, 6mm diameter
936200	Rod 200mm long, 6mm diameter

CLAMPS*

Code	Description
93010	Large Clamp
93110	Medium Clamp
93310	Small Clamp
93020	Multi-screw Clamp
93030	Large-Medium Transition Clamp
93120	Medium Multi-screw Clamp
99-93040	Large Double Multiscrew Clamp
99-93140	Medium Double Multiscrew Clamp

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The Orthofix Galaxy Fixation System components not listed above have not been tested for heating, migration, or image artifact in the MR environment, and their safety is unknown. Scanning a patient carrying a frame that includes these components may result in patient injury.

* Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Orthofix representative if you have questions about the availability of Orthofix products in your area.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

Symbol	Description	
 	Consult instructions for Use	
	Single use. Do not reuse	
	STERILE. Sterilised by irradiation	
	NON STERILE	
REF	LOT Catalogue number	
	Expiry date (year-month-day)	
CE	CE 0123 CE marking in conformity to applicable European Directives/Regulations	
	 Date of manufacture	
	Do not use if package is opened or damaged	
	Symbol for MR conditional. This means an item has been demonstrated to pose no known hazards in a specified MRI environment with specified conditions of use.	
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.	

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